



Medical

Pall Lipipor™ TNA Filter Set



Air-eliminating filter set for lipid-containing nutritional IV administration

Features

- ▶ Pall Lipipor 1.2 µm nylon membrane
- ▶ Smooth housing design
- ▶ Sterile, non-pyrogenic fluid pathway
- ▶ Phthalate-free tubing extension

Benefits

- ▶ Removes particles and microorganisms
- ▶ Retains oversized lipid droplets
- ▶ Eliminates air (independent of filter position)
- ▶ Easy to prime and use
- ▶ Increases patient comfort